

Aseptic Technique



A GMP/GTP Training Module




Aseptic Technique



- The GMP Facility manufactures products for clinical use
- These products must meet a number of requirements, one of which is that they are **sterile**
- In order to produce sterile components it is necessary to follow **aseptic technique**

Aseptic Technique




- **Use of aseptic technique protects the product and the operator from exposure to infectious agents and contaminants**
- **The technique is a way of working that provides this protection**
- **The following presentation provides basic information on aseptic technique**

Aseptic Technique



- **The idea behind aseptic technique is to PREVENT contamination**
- **It is not possible to “clean up” an already contaminated product by using aseptic technique**
- **In order to maintain sterility, it is necessary to start with a sterile product & handle it in an aseptic manner**

Before Starting



- **All reagents and materials to be used should be in-date, barcoded, and should have a Certificate of Analysis on file**
- **All equipment to be used should have been cleaned according to SOPs and that cleaning should be documented**
- **All equipment, including pipettors, should have been calibrated**
- **Waste traps should contain some bleach to act as a decontaminant**

Before Starting



- **Collect all of the reagents and materials that you will need onto a workspace close to the biological safety cabinet**
- **Make sure that you have a biohazard trash container, sharps container and waste trap available – we recommend that waste traps be located on the floor OUTSIDE the BSC**
- **Make sure that you have disinfectant spray bottles containing 3M Quat and 70% ethanol available**

Before Starting

- **Make sure that you are wearing the required personal protective equipment**
 - **Blue labcoat**
 - **Gloves**
 - **Eye protection if necessary**
 - **Sleeve covers available (if you will work with more than one product)**




Before Starting

- Page QA/QC if you require environmental monitoring during production
- Use alcohol foam to clean gloves before using the BSC



Before Starting

- Wipe down the working surface in the BSC with 3M Quat or 70% ethanol and leave to dry
- Disinfect any equipment e.g. racks, before placing them into the BSC



Starting Work

- Do not clutter up the working area in the BSC – bring in only what is necessary
- Be careful not to block the air vents at the back and front of the BSC



Starting Work

- When you start working, be sure that you work well inside the BSC and not directly above the front air vents
- Keep your nose and mouth away from the window opening





Working Tips




- **Tighten the caps on culture flasks before moving them to the BSC**
- **Wipe down culture flasks with spray disinfectant before bringing them into the BSC**
- **Decontaminate your gloves using alcohol foam frequently during the culture procedure**

Working Tips

- Be careful not to touch the necks of culture flasks when removing the caps
- Be careful not to touch the inside of the cap
- Place the cap either face up, or face down on an alcohol wipe



Working Tips



- **Handle cells gently when pipetting and resuspending**
- **Do not over-expose cells to enzymes when dispersing a monolayer into a suspension**
- **Avoid extreme temperature changes**

Working Tips

- Resuspend in medium that can maintain the correct pH of the culture
- Bicarbonate-buffered medium requires a 5% CO₂ atmosphere to maintain its pH
- Most media contain pH indicators and are a reddish orange at the correct pH



Working Tips



- **A culture that becomes cloudy rapidly after incubating is probably contaminated**
- **This can be checked by sterility cultures and microscopic examination**
- **QA/QC recommends that all contaminated cultures should be sterility tested to identify the contaminating organism**
- **Cultures that contain many cells usually turn yellow, but the supernatant medium remains clear rather than cloudy, as with contaminated cultures**

Working Tips



- Cultures that contain white or gray threads, or cotton wool-like floating particles usually are contaminated with fungus or mold
- This can be checked by sterility cultures and microscopic examination
- QA/QC recommends that all contaminated cultures should be sterility tested to identify the contaminating organism
- Fungus and mold is difficult to cure effectively

Working Tips



- **Contaminated cultures should normally be discarded after a sterility culture has been submitted**
- **If it is essential to maintain a contaminated culture, it should be moved to a separate quarantine incubator – see QA/QC**
- **The incubator in which it was grown should also be cleaned, including wiping down the exterior of any flasks that were sharing the incubator**

Working Tips



- **When returning a flask from the BSC to the incubator, tighten the cap first**
- **When the flask is in the incubator, the cap may then be loosened to allow for gas exchange**
- **Avoid medium entering the neck of the flask during transport**
- **Any medium that spills in the incubator should be sprayed with disinfectant and immediately wiped up**

Working Tips



- **QA/QC recommends that medium that is prepared/supplemented for cell culture should be sterility tested before use whenever possible**
- **Random bottles from new lots should be submitted to QA/QC for testing as a check on the manufacturer**
- **Remember to assign an expiration date to newly supplemented medium, and to indicate the nature and amount of the supplement(s)**
- **Initial and date the changes to the label.**

Serological Pipets

When working with serological pipets

- Transfer pipet to the BSC before opening
- Open package by grasping each side of the closure and peeling back the wrapper from the wide end
- Do not touch the unwrapped area of the pipet
- Plug the unwrapped end firmly into the Pipet Aid device



Serological Pipets

When working with serological pipets

- Remove the wrapper by pulling it from the tip end of the pipet, taking care not to touch the pipet



Serological Pipets

- **Be careful not to let the pipet touch items inside the BSC**
- **When pipetting liquids, insert the pipet tip beneath the surface of the liquid, taking care not to let the pipet barrel touch the neck of the container**




Serological Pipets



- Draw up the liquid into the barrel, do not blow air thru the liquid
- Remove the pipet from the container, again avoiding touching the barrel against the container neck

Serological Pipets



- **Take care not to let liquid drip from the pipet when transferring**
- **Insert the pipet into the receiving vessel without touching the barrel to the neck**
- **Gently expel the liquid into the vessel. Do not blow bubbles. Expel the liquid against the vessel wall if necessary**
- **The same pipet may be used to resuspend cells in the liquid, but must then be discarded**

Serological Pipets



- Remove the pipet from the Pipet Aid by grasping the barrel firmly and pulling
- Discard the pipet into the SHARPS container
- Disinfect gloves with alcohol foam
- If any liquid has dripped onto surfaces, clean it up with an alcohol swab, or disinfectant spray and gauze.



Micro-Pipetting



- When using Eppendorf or similar pipets, individually wrapped tips should be handled like serological pipets
- Racked pipet tips can be used, but care must be taken to keep the tip box covered when tips are not being removed
- **If you are in doubt about the sterility of a rack of pipet tips, discard them, and use a new rack**
- Wipe down pipet barrels after use using an alcohol wipe

Micro-Pipetting




- **The use of barrier tips is recommended for most applications**
- **If liquid enters the pipet barrel, stop using the pipetter until it can be properly cleaned and disinfected**
- **When using micro-pipettors for cell sampling, be sure to mix the cell suspension thoroughly**

Universal Precautions




- **All cell products should be regarded as potentially infectious**
- **Follow Universal Precautions when working with these products**
- **You will receive annual training on Universal Precautions**

General Precautions



- **Avoid creating aerosols when pipetting or transferring liquids**
- **Handle cells gently – do not shake or pipet suspensions vigorously**
- **If you feel that you may have touched a sterile item – discard it and use a new one**
- **Clean up any spills immediately after decontaminating with disinfectant**

General Precautions



- **Make sure caps are on tight before centrifuging cell suspensions**
- **If possible use aerosol containment or wrap parafilm around the caps**
- **Avoid spinning cells too hard or for too long**
- **Pellet resuspension is best done gently, using small volumes of medium, followed by dilution of the suspension**

Waste Reagents

- Vacuum aspiration is the easiest way to remove supernatants [SN]
- The SN should be aspirated into a waste container containing undiluted bleach
- If the container is not full at the end of the day, be sure to add more bleach – the liquid in an active waste container should normally be a clear yellow



Waste



- **DO NOT empty full liquid waste containers**
- **The container should be disconnected from the vacuum apparatus and discarded into the biohazardous trash**
- **Discard old culture flasks into the biohazardous trash after tightening the caps**
- **All pipets should be discarded into the SHARPS containers**
- **DO NOT overfill waste containers**